

# 100 bp DNA Ladder Plus (100bp-3kb) Ready-to-Use

**Description:** A DNA ladder composed of 12 fragments ranging from 100bp - 3kb. The intensity of the 500bp, 1.5kb bands have been increased to serve as an internal reference.

**Storage:** -20° C for up to 24 months, or 4° C for up to 3 months.


**Recommend loading:** 5 µl (or 1 µl/mm width of gel)

**Concentration:** 108 µg/ml

SKU	Size	Use
DNA909.50	50 µg	100 applications
DNA909.550	5 x 50 µg	500 applications

**DNA Mass**  
(ng/5 µl)

**Base Pairs**



1.5% TAE Agarose gel